



MRidium™*
Disposable Products for MR Infusion Pump

The delivery of critical medications during MRI procedures increasingly requires simultaneous infusion of multiple medications, that there be minimal connection and disconnection from the patient's IV site, and that patient set-up be fast and efficient. Covidien provides the MRidium™* MR infusion system bypass infusion set, standard infusion set and syringe set to aid this process.

Order Information

Description	Qty	Unit	Size	Order #	11 Digit NDC
Tubing, infusion, bypass	50	set	~215 cm (85")	MI1055	
Tubing, infusion, standard	50	set	~300 cm (118")	MI1056	
Tubing, syringe	50	set	~216 cm (85")	MI1057	

Features and Benefits

Feature	Benefit
Unique bypass infusion set	Provides a quick and easy means of switching from a hospital pump to the MRidium™* MR infusion system without disconnecting the IV from the patient. The bypass set minimizes the need for long-tubing workarounds that may affect the delivery of patient medications.
Range of tubing set priming volumes	Volume capabilities include ~ 4 mL, ~12 mL and ~18 mL to meet a variety of needs

Frequently Asked Questions

Question	Answer
Can I use other brands of infusion sets with MRidium™* MR infusion system?	Every pump uses its own proprietary tubing. The MRidium™* MR infusion system utilizes only MRidium™* infusion sets. The bypass set will attach to your hospital tubing at the upper and lower Y sites, making a circuit around the original hospital tubing so it does not need to be discarded or disconnected from the patient.

*MRidium is a trademark of Iradimed Corporation and used under license.
