



Optivantage™

Disposable Products for Dual-Head CT Contrast Delivery System

The Optivantage™ dual-head CT syringe is 200mL, front-loading, sterile and disposable. The semi-transparent syringe is designed to aid visual inspection for small air bubbles that may be present prior to completing the injection. We offer CT multipacks in a variety of combinations to fit your individual requirements, with Optivantage™ Y-tubing in dual, single and no check valve options.

Order Information

Description	Qty	Unit	Size	Order #	11 Digit NDC
Multipack, with syringe and low pressure coiled tubing	50	each	200 mL, 60"	800099	
Multipack, with syringe and Y-tubing, dual check valve	50	each	200 mL, 60"	844021	
Multipack, with syringe and Y-tubing, no check valve	50	each	200 mL, 60"	844020	
Multipack, with syringe and Y-tubing, single check valve	50	each	200 mL, 60"	844015	
Multipack, with two syringes and 60" Y-tubing, dual check valve	20	each	200 mL, 60"	844023	
Multipack, with two syringes and Y-tubing, no check valve	20	each	200 mL, 60"	844022	
Multipack, with two syringes and Y-tubing, single check valve	20	each	200 mL, 60"	844016	
Syringe	50	each	200 mL	800096	
Transfer set, CT	50	each		810551	
Tubing, coiled extension, low pressure	100	each	60"	601195	
Tubing, straight extension, low pressure	100	each	48"	601277	
Y-tubing, coiled, no check valve	50	each	60"	844010	
Y-tubing, dual check valves	50	each	60"	844012	
Y-tubing, single check valve	50	each	60"	844011	

Features and Benefits

Feature	Benefit
Semi-transparent design	Optivantage™ syringes are designed to aid visual inspection for air bubbles that may be present prior to completing the injection, which may reduce the chance of injecting air and may aid patient safety.
Package Combinations	Optivantage™ disposable syringes and Y-tubing options are offered in cases of single, individually packaged pieces. If preferred, the multipack combination provides the components needed to perform one procedure in a single tray.
Transfer Set	The transfer set provides quick and easy saline filling, aiding efficiency for dual head procedures.

**Multiple Y-Tubing
Combinations**

Our Y-tubing is available with no check valve, single check valve and dual check valves for maximum flexibility to meet your individual needs and requirements.

Frequently Asked Questions

Question

Answer

**What does the Optivantage™
CT multipack contain?**

Each contains 60" coiled Y-tubing, a Handi-Fil™ fill straw, and one or two 200mL empty syringe (s). Customers using Ultraject™ prefilled syringes would need multipacks containing only one 200mL syringe.

**Why single vs. dual check
valves on the 60" coiled Y-
tubing?**

Dual check valves offer the maximum protection against the mixing of contrast and saline in the syringes. A single check valve offers mixing protection for the contrast syringe while still allowing the clinician to "draw back" if desired.

Instructions For Use

Instructions for use may be found in the Optivantage™ Operator's Manual.
