



CT 9000™ ADV

Disposable Products for CT Contrast Delivery System

The CT 9000™ ADV syringe is 200mL, front-loading, sterile and disposable. The semi-transparent syringe is designed to aid visual inspection for small air bubbles that may be present prior to completing the injection. We offer 60" coiled low pressure extension tubing and 48" straight low pressure extension tubing. The syringe and 60" coiled tubing are also combined in the CT 9000™ ADV multipack, providing convenient options to meet your needs.

Order Information

Description	Qty	Unit	Size	Order #	11 Digit NDC
Multipack, with syringe and low pressure coiled tubing	50	each	200 mL, 60"	800099	
Syringe	50	each	200 mL	800096	
Transfer set, CT	50	each		810551	
Tubing, coiled extension, low pressure	100	each	60"	601195	
Tubing, straight extension, low pressure	100	each	48"	601277	
Y-tubing adaptor for connecting Optistat™ contrast delivery system to CT 9000™ ADV contrast delivery system	50	each	6" (each adaptor leg)	810555	

Features and Benefits

Feature	Benefit
CT 9000™ ADV multipack	Convenient multipacks provide the products needed for your procedure type in one simple tray
Semi-transparent	CT 9000™ ADV syringes are designed to aid visual inspection for air bubbles that may be present prior to completing the injection, which may reduce the chance of injecting air and may aid patient safety.

Frequently Asked Questions

Question	Answer
What is included with CT 9000™ ADV multipack?	The CT 9000™ ADV multipack includes one 200 mL syringe, one Handi-Fil™ fill straw and one 60" low pressure coiled tubing.

Instructions for Use

Instructions for use may be found in the CT 9000™ ADV Operator's Manual.
