



MRidium™*
Accessories for MR Infusion System

The MRidium™* MR infusion system is the first and only non-magnetic infusion pump for the MR suite. Optional accessories include the wireless remote control, giving clinicians the freedom to operate the system from the control room, and the MRidium™* SideCar™* module for the second pump channel, allowing critically ill patients to remain on one or two IV medications, before, during, and after a MRI procedure.

Order Information

Description	Qty	Unit	Size	Order #	11 Digit NDC
Pole, IV, non-magnetic	1	each		MI1119	
Pump, non-magnetic infusion, Sidecar™ for channel B side	1	each		MI3851	
Remote control for non-magnetic infusion pump	1	each		MI3855	

Features and Benefits

Feature	Benefit
Second-channel SideCar™* module	The second-channel SideCar™* module allows critically ill patients to remain on one or two IV medications, before, during, and after a MRI procedure. Continuous infusions facilitate uniform drug delivery throughout the procedure. Each channel can deliver medications via syringe, bag, or bottle, based on clinician needs.
Wireless remote control	The optional wireless remote control gives clinicians the freedom to operate the system from the control room – outside the MR suite. The remote control can be used to monitor and control all dosing parameters, including titration of patient meds. Please note: The remote control cannot be used in the MR suite.

Frequently Asked Questions

Question	Answer
Can the remote control be used in the MR suite?	No
Is the IV pole non-magnetic?	Yes

*MRidium and SideCar are trademarks of Iradimed Corporation and used under license.
